

Cylinder holder

Pack contents and wall mounting instructions.

1. Choose suitable screws for the wall surface that the cylinder holder is to be mounted.
2. Ensure the wall surface that the holder is to be mounted is clean and flat:
 - if existing surface is not suitable, consider mounting onto a stable mounting board.
3. Do not place in outdoor and high moisture environments.
4. Ensure that the back plate covering the springs is in place before securing the holder to the wall.
5. Screws should be no larger than 6mm in diameter:
 - choose screws with pan heads no larger than 12mm diameter;
 - use the included template for accurate screw positions.
6. Please follow manufacturer's instructions when determining appropriate screw depth.
7. Allow 22mm max for mould depth when choosing screw length.
8. A total of four screws are required.
9. Mount a minimum of 200mm from the floor.
10. Successive units should be mounted no closer than 185mm spacing from the bracket centre line.
11. The strap has been supplied to provide additional cylinder security if required.
 - shoulder screws, for use with the strap, have 'Patlok' coating on the threads to allow semi-permanent fixing;
 - if the strap is not to be used, the enclosed paddle plugs can be pushed into the paddle inserts.



INALO setup

To find out more or place your order today please call
1300 363 109 or visit www.boc.com.au

Australia: BOC Limited ABN 95 000 029 729
BOC Healthcare, 10 Julius Ave, North Ryde NSW 2113
T: 1300 363 109 F: 1300 363 438 E: hospital.care@boc.com

Details given in this document are believed to be correct at the time of printing. While proper care has been taken in the preparation, no liability for injury or damage resulting from its use can be accepted. BOC is a trading name of BOC Limited, a member of The Linde Group. © BOC Limited 2012. Reproduction without permission is strictly prohibited.